**MEDICAL**  
TEMPLATES

# PUNCTURE CUBE

**MANUALE D'USO**  
USER MANUAL

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
<b>USER MANUAL.....</b>	<b>18</b>
<b>1 - INTRODUCTION.....</b>	<b>18</b>
<b>1.1 - INTENDED USE.....</b>	<b>18</b>
<b>1.2 - INTENDED PURPOSE.....</b>	<b>19</b>
<b>1.3 - PRODUCT OVERVIEW.....</b>	<b>20</b>
<b>2 - IMPORTANT INSTRUCTIONS.....</b>	<b>22</b>
<b>2.1 - MEANING OF THE SYMBOLS.....</b>	<b>22</b>
<b>2.2 - IMPORTANT SAFETY INSTRUCTIONS.....</b>	<b>22</b>
<b>2.3 - SELECTING THE PUNCTURE INSTRUMENT.....</b>	<b>23</b>
<b>2.4 - PATIENT SAFETY INSTRUCTIONS.....</b>	<b>23</b>
<b>3 - APPLICATION.....</b>	<b>24</b>
<b>3.1 - PREPARATION.....</b>	<b>24</b>
<b>3.2 - PLACEMENT OF THE PUNCTURE CUBE.....</b>	<b>25</b>
<b>3.3 - IDENTIFICATION OF THE COORDINATES.....</b>	<b>25</b>
<b>3.4 - PLACING THE PUNCTURE NEEDLE.....</b>	<b>26</b>
<b>3.5 - COLLAPSING THE PUNCTURE CUBE.....</b>	<b>27</b>
<b>3.6 - PERFORMING THE TREATMENT.....</b>	<b>28</b>
<b>3.7 - REMOVING THE PUNCTURE CUBE.....</b>	<b>28</b>

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


# SHORT GUIDE

An aid for needle alignment during CT-guided interventions, for example infiltrations, injections, biopsies, or placing access points for drainages. Can be used anywhere on the body, on adults and children, provided that it is properly adhered to the patient. The device must be used in combination with a CT and associated software, by a physician and trained personnel.

## 1. Place the Cube

Place the PUNCTURE CUBE over the target region of the intervention. Attach to the patient, ensuring all four adhesive tabs are firmly adhered to the skin.

## 2. Scan

Perform a standard planning scan, being sure to capture at least 50% of the PUNCTURE CUBE. Load the scan in the planning software.

## 3. Detect the Cube

Use the Recognize Cube function and activate the virtual needle in the planning software.

## 4. Plan the Intervention

Use the virtual needle in the software to plan a trajectory. Note the final coordinates on the upper and lower plates of the PUNCTURE CUBE as suggested by the software. Be sure to note the corner as well.

## 5. Perform the Puncture

Insert the needle through the star-shaped holes on the upper and lower plates of the PUNCTURE CUBE. Lean the needle into the corner as indicated by the software. Check the trajectory with control scans. Remove the clips and lift the top plate of the PUNCTURE CUBE in order to collapse the Cube if necessary.

This short guide does not replace the full user manual. The user must be familiar with the procedure for safe use of the product.

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# USER MANUAL

Thank you for purchasing our device. The PUNCTURE CUBE is characterized by its efficiency and ease of use. Please read the instruction manual carefully to ensure that you get the most out of the PUNCTURE CUBE. Please contact Medical Templates or your local medical device adviser if you have any additional questions about this device or how to use it.

## MANUFACTURER INFORMATION

MEDICAL TEMPLATES AG  
Brunnenwiese 10  
8132 Egg bei Zürich  
Switzerland

[www.medicaltemplates.ch](http://www.medicaltemplates.ch)  
[info@medicaltemplates.ch](mailto:info@medicaltemplates.ch)

## 1 - INTRODUCTION

This instruction manual forms an integral part of the device. It describes its intended use and must be kept near the device. Following the PUNCTURE CUBE instruction manual is essential for handling the device properly and making sure it works correctly. It ensures the safety of the patient and user alike.

### 1.1 - INTENDED USE

The PUNCTURE CUBE is used for guiding needles and instruments during Computed Tomography (CT)-guided punctures, such as infiltrations, diagnostic and therapeutic biopsies, or for placing access points for drainages or ablations. The PUNCTURE CUBE can be used anywhere on the body, on adults and children, provided that the PUNCTURE CUBE is resting on all four feet and is properly adhered to the patient. The device must be used in combination with a CT and associated software by a physician. The PUNCTURE CUBE is a single-use sterile device.

## **1.2 - INTENDED PURPOSE**

### **Medical Indication**

The PUNCTURE CUBE is used for guiding needles and instruments during CT-guided punctures, such as infiltrations, diagnostic and therapeutic biopsies, or for placing access points for drainages or ablations. The following cannula and needle diameters are permitted for use with the PUNCTURE CUBE: G18 - G22.

### **Patient Profile**

Adults and children from 3 years and up. It can be used anywhere on the body, provided that the PUNCTURE CUBE is resting on all four feet and is properly adhered to the patient.

### **User profile**

The PUNCTURE CUBE can be unpacked and attached to the patient by a radiographer. A trained physician must perform all other functions.

### 1.3 - PRODUCT OVERVIEW


Fig. 1:  
PUNCTURE CUBE in its  
packaging  
(Top and bottom view)



Fig. 2:  
PUNCTURE CUBE  
unpacked


Fig. 3


Fig. 4

The PUNCTURE CUBE consists of two plastic components: an upper template (Fig. 3) and a lower template (Fig. 4). Both templates consist of a grid of star-shaped openings. The grids are divided into quadrants containing twelve openings each. The markings on the sides of the grid clearly identify each quadrant.


Fig. 5


Fig. 6

Four adhesive tabs, covered with removable protective film, are attached to the lower template. The tabs are used to attach the PUNCTURE CUBE to the patient. The upper and lower template are joined together into a cube. The upper template is held in place with two clips to avoid unintended disassembly (Fig. 5 and Fig. 6).

The PUNCTURE CUBE is affixed to the body fully assembled.


Fig. 7


Fig. 8

The star-shaped openings on the lower and upper templates have four corners which serve as guides for the needle (Fig. 7). When using the PUNCTURE CUBE, the needle should lean into the indicated corner of the given coordinate on both the upper and lower templates (Fig. 8).

## 2 - IMPORTANT INSTRUCTIONS

### 2.1 - MEANING OF THE SYMBOLS

	Manufacturer
	Lot number
	Reference number
	Date of manufacture
	Use by date
 0297	This device complies with the essential requirements set out in Directive 93/42/EEC concerning medical devices
	Consult instructions for use
	Sterilized using ethylene oxide
	Keep away from sunlight
	Keep dry
	Storage temperature: +10°C to +27°C
	Do not use if package is damaged
	Do not re-use
	Authorized representative in the European Community
	Importer in the European Community
	Medical device

### 2.2 - IMPORTANT SAFETY INSTRUCTIONS

Important text passages that require special attention are marked with a warning triangle.

The PUNCTURE CUBE is packaged and delivered in assembled condition.

The PUNCTURE CUBE should arrive fully assembled. If the templates shift during transportation, the following steps must be followed for manual assembly:

- First, fold the two long sides of the lower template up at a right angle and hold them firmly.
- Fold the two short sides upwards. Attach the sides together by inserting the tabs of the long sides into the slots of the short sides.
- Place the upper template on top. When aligning, ensure that the arrows on both templates point in the same direction.
- Gently press the top template so that the tabs of the side panels fit into the slots of the the template. Snap the provided clips onto the frame to ensure the PUNCTURE CUBE is firmly assembled.

The PUNCTURE CUBE is a sterile, single-use device. The PUNCTURE CUBE may only be handled when wearing sterile disposable gloves.

## **2.3 – SELECTING THE PUNCTURE INSTRUMENT**

In order to ensure proper functioning of the PUNCTURE CUBE, the following sizes of cannula and needle are permissible for use: G18, G19, G20, G21, G22.

When selecting the needle, it must be noted that the height of the PUNCTURE CUBE in the collapsed state is 12 mm. This distance must be added to the total penetration depth when determining the required length of the needle. Also note that when the angle of incidence is extremely diagonal, the available needle length is shortened considerably.

## **2.4 – PATIENT SAFETY INSTRUCTIONS**

The PUNCTURE CUBE is affixed to the patient with adhesive tabs. The adhesive, Acrylate, is approved for medical use, but has a risk of causing slight skin irritation.

Always perform a control scan to ensure correct positioning of the needle, as the target area can move after the initial scan due to breathing.

Always ensure the needle is inserted through the holes and corners indicated by the software. Using the wrong coordinates can lead to a missed target.

When inserting a needle or thin instrument through the PUNCTURE CUBE, there is a small risk that part of the PUNCTURE CUBE could be scraped by the needle and inserted into the patient. Place the needle carefully to avoid this.

## **3 – APPLICATION**

### **3.1 – PREPARATION**

- ① Perform a topogram to identify the scan area required for the puncture and mark the spot where the PUNCTURE CUBE will be attached.
- ② Prepare the site where the PUNCTURE CUBE will be attached: Clean and disinfect the site as well as the surrounding area. Ideally, place sterile drapes around the PUNCTURE CUBE, so that the sides fall within the sterilized area when collapsed.
- ③ Carefully peel off the protective film from the packaging and place it aside, as it will be needed again at a later stage.
- ④ Wearing sterile gloves, remove the PUNCTURE CUBE from the package and pull the protective film covers off the adhesive tabs.

### 3.2 – PLACEMENT OF THE PUNCTURE CUBE

- ① Position the Puncture Cube over the patient and orient the it correctly.

The top template's arrow must point towards the cranium. The PUNCTURE CUBE can be mounted transversally or longitudinally.

- ② Place the PUNCTURE CUBE on the patient's body in the desired location and press the PUNCTURE CUBE gently.
- ③ Press the adhesive tabs firmly against the skin to ensure the PUNCTURE CUBE is securely attached. All four feet of the PUNCTURE CUBE must be firmly adhered to the patient's skin. Insufficient fixation of the feet can lead to the PUNCTURE CUBE shifting after the planning scan and an inaccurate puncture.

Once the adhesive tabs have been pressed on, repositioning is no longer possible!

### 3.3 – IDENTIFICATION OF THE COORDINATES

- ① Perform a spiral CT scan of the PUNCTURE CUBE.

At least half of the PUNCTURE CUBE must be captured. The PUNCTURE CUBE may not be covered or taped over at any point during the CT scan. Only then is the identification of the PUNCTURE CUBE by the planning software ensured.

- ② Load the recorded images in compatible planning software to register the PUNCTURE CUBE and plan the puncture. Please refer to the planning software instruction manual for this purpose. When using planning software from MEDICAL TEMPLATES / Synedra, the user manual is included.

- ③ Note down the coordinates provided by the planning software on the form provided for this purpose, which can be found on the sealing film.

Ensure that the PUNCTURE CUBE's position has not changed between the planning scan and the puncture.

### 3.4 – PLACING THE PUNCTURE NEEDLE

- ① Locate the coordinates of the star-shaped hole specified by the software on the top template of the PUNCTURE CUBE with the help of guidelines A-D and 1-4. Position the puncture needle in the selected star-shaped hole, making sure the needle is also resting in the corner of the hole specified in the software.
- ② Position the needle in the bottom template using the same principle. Here, the guidelines contain the letters E-G and 5-7.

The markers next to the grid in the form of triangle, semicircle, circle, and square are only used for identification by the software and are not used for puncture.

Local anesthesia can now be administered at the injection site. To do this, place an anesthetic needle vertically through any two holes in the PUNCTURE CUBE to reach the injection site.

- ③ Ensure the needle is inserted through the holes and corners indicated by the software. Perform the planned puncture with the positioned needle.

Do not apply pressure to the PUNCTURE CUBE while inserting the needle, to ensure you do not change the position of the PUNCTURE CUBE. This is the only way to ensure that the needle position calculated by the software can be maintained.

While inserting the needle, the pressure on the needle must not be so high that it causes the needle to bend. Only apply pressure in the direction of the needle tip.

When using needles with beveled tips, make sure that the needle follows a straight trajectory (e.g. rotate the beveled needle while piercing).

Finish the first puncture at least 2 cm before the proximal end of the needle meets the top template. Then the top template can be lifted off and then placed on the lower template if needed.

- ④ Perform a control scan to check the needle's position. Repeated control scans may be necessary to insert the needle safely.

### 3.5 – COLLAPSING THE PUNCTURE CUBE

If necessary, the PUNCTURE CUBE can be collapsed. This allows the needle to reach the required penetration depth without restriction by the PUNCTURE CUBE.

If doing so is not necessary, continue to step 3.6 - Performing the treatment.

To unfold the PUNCTURE CUBE:

- Remove the clips on the top template while firmly holding the PUNCTURE CUBE on two sides.
- Carefully lift the upper template until the tabs have loosened from the slots. The side sections will fold out to the side automatically.
- Slide the top template down along the needle and place it onto the lower template.

Do not pull the upper template upwards, otherwise the needle will be retracted

Continue to guide the needle down to the target. Additional control scans may be necessary.

### **3.6 – PERFORMING THE TREATMENT**

Perform the necessary treatment (e.g. injection or aspiration).

### **3.7 – REMOVING THE PUNCTURE CUBE**

Carefully remove the needle from the patient and PUNCTURE CUBE after treatment.

Detach the PUNCTURE CUBE from the skin of the patient and dispose of the PUNCTURE CUBE.

Reuse of the PUNCTURE CUBE is not permitted.

In case of further questions, please contact our support team:  
[support@medicaltemplates.ch](mailto:support@medicaltemplates.ch)

MEDICAL TEMPLATES AG  
Brunnenwiese 10  
8132 Egg, Zurigo / Egg, Zürich  
Svizzera / Switzerland

[www.medicaltemplates.ch](http://www.medicaltemplates.ch)  
[info@medicaltemplates.ch](mailto:info@medicaltemplates.ch)

Numero di riferimento / Reference Nr.: mt-pc-ifu-v56-it  
Aggiornamento / Update: 02/03/2023  
[www.medicaltemplates.ch/ifu](http://www.medicaltemplates.ch/ifu)

MedEnvoy Global B.V.  
Prinses Margrietplantsoen 33 - Suite 123  
2595 AM The Hague  
Paesi Bassi / The Netherlands